

MYO PC SOFTWARE



INSTALLATION

- CD starts automatically
- Choose the Installation option in the main menu



MINIMUM CONFIGURATION FOR PC

- Operating System: Windows 7, 8 & 10
- CPU: Intel Core i5 - 1,60 GHz
- System Memory (RAM): 2048MB
- Hard disc: 1GB free space
- Screen: minimum 1280 x 800
- USB2



GENERAL SAFETY WARNING

- Only qualified people who are trained in the application may use this software.
- Before installing the software, read the instructions in the user manual available on this CD-ROM.
- Make sure that you understand the functionality of the Myo 200 before using it together with the MYO PC software.
- Read the user instructions of the Myo 200.
- Install this software on a personal computer which function is uncompromised.
- Do not use the PC for other purposes during the treatment.
- No other accessories should be connected during treatment that could compromise the functionality of the software.
- If there is any doubt on the functionality of the software, contact the manufacturer.
- Always keep the user instructions with the equipment.

ELPC01 2-GAM 02-V4.1



GymnaUniphy NV | Pasweg 6a | B-3740 Bilzen | Belgium
Tel +32(0)89 510.510 | Fax +32(0)89 510.511
info@gymna.com | www.gymna.com

CE 0344

gymna